



## Disposable Auto-Feed Adult/Pediatric Humidification Chamber VHC10

Ref	VHC10
Inlet / Outlet(mm)	22M / 22M
Compressible Volume(Empty)	300 ml
Compressible Volume(Full)	140 ml
Maximum Water Capacity	160 ml
Maximum Peak Flow	90 l/min
Heat Conduction Diameter	12.2cm
Resistance (cmH2O)	(60L/Min), < 1cmH2O
Compliance(Empty)	≤ 0.41ml/cmH2O
Package	12pcs/carton



## Disposable Auto-Feed Neonatal Humidification Chamber VHC20

Ref	VHC20
Inlet / Outlet(mm)	22M / 22M
Compressible Volume(Empty)	430 ml
Compressible Volume(Full)	195 ml
Maximum Water Capacity	235 ml
Maximum Peak Flow	90 l/min
Heat Conduction Diameter	12.2cm
Resistance (cmH2O)	(60L/Min), < 1cmH2O
Compliance(Empty)	< 0.7ml/ cmH2O
Package	12pcs/carton

## Reusable Humidification Chamber VHC25



Ref	VHC25
Inlet / Outlet(mm)	22M / 22M
Compressible Volume(Empty)	640 ml
Compressible Volume(Full)	140 ml
Maximum Water Capacity	500 ml
Maximum Peak Flow	90 l/min
Heat Conduction Diameter	12.2cm
Resistance (cmH2O)	(60L/Min), <1cmH2O
Compliance(Empty)	< 0.85ml/ cmH2O
Package	12pcs/carton