



## HEPA

FSC No.	504-010
Part No.	51050947
Connector(mm)	22F/15M-22M/15F
Use	Disposable
Media Type	"Hydrophobic mechanical glass-fibre"
NaCl Filter Efficiency	99.97%
Pressure Drop	2.5cmH2O@30L/min
Compressible Volume	65ml
Weight	48g
Package	100pcs/carton



## Low Resistance bacterial / Viral respiratory filter BSF 103-104

FSC No.	510-014	510-014
Part No.	51008900	51009000
Compressible Volume	65ml	65ml
Pressure Drop at 30LPM	1cmH2O	1cmH2O
Bacterial Filtration Efficiency	99.998+%	99.9993+%
Viral Filtration Efficiency	99.995+%	99.996+%
Casing Material	k-resin	k-resin
Weight	27g	27g
Package	100pcs/carton	100pcs/carton

## High performance bacterial / Viral respiratory filter BSF 201-204



FSC No.	504-011	504-011	504-011	504-011
Part No.	50421500	50421700	50452200	50453700
Pressure Drop at 30LPM	1.5cmH2O	1.5cmH2O	1.5cmH2O	1.5cmH2O
Bacterial Filtration Efficiency	99.99994+%	99.99994+%	99.99994+%	99.99994+%
Viral Filtration Efficiency	99.997+%	99.997+%	99.9999+%	99.9999+%
Gas Sampling port	yes	no	yes	No
Weight	25g	25g	25g	25g
package	100pcs/carton	100pcs/carton	100pcs/carton	100pcs/carton

## HMEF Series



### HMEF 1000 Straight

FSC No.	504-011
Part No.	50497100
Connector	22F/15M-22M/15F
Dead Space	35ml
Tidal Volume Range	150~1200ml
Pressure Drop	0.35kPa@30L/min
Moisture Output	9.5mg/L@Vt150~1200ml
Bacterial Filtration Efficiency	>99.9999%
Viral Filtration Efficiency	>99.99%
Weight	30g
Package	100pcs/carton



### HMEF Midi Straight

FSC No.	504-011
Part No.	50400216
Connector	22F/15M-22M/15F
Dead Space	20ml
Tidal Volume Range	60~500ml
Pressure Drop	0.35kPa@30L/min
Moisture Output	9.5mg/L@Vt60~500ml
Bacterial Filtration Efficiency	>99.99%
Viral Filtration Efficiency	>99.99%
Weight	20g
Package	100pcs/carton



### HMEF 1000 Angle

FSC No.	504-011
Part No.	50400218
Connector	22F/15M-22M/15F
Dead Space	33.5ml
Tidal Volume Range	150~1200ml
Pressure Drop	0.35kPa@30L/min
Moisture Output	9.5mg/L@Vt150~1200ml
Bacterial Filtration Efficiency	>99.9999%
Viral Filtration Efficiency	>99.99%
Weight	30g
Package	100pcs/carton



### HMEF Midi Angle

FSC No.	504-011
Part No.	50400216
Connector	15M-22M/15F
Dead Space	24ml
Tidal Volume Range	60~500ml
Pressure Drop	0.35kPa@30L/min
Moisture Output	9.5mg/L@Vt60~500ml
Bacterial Filtration Efficiency	>99.99%
Viral Filtration Efficiency	>99.99%
Weight	20g
Package	100pcs/carton

## HME Trach



### HME Trach Basic

FSC No.	510-019
Part No.	5150901
Connector	15F
Dead Space	14ml
Pressure Drop	0.03KPa@30L/min
Moisture Output	37mgH <sub>2</sub> O/L
Weight	15g
Package	200pcs/carton



### HME Trach T

FSC No.	510-019
Part No.	5150902
Connector	15F
Dead Space	8ml
Pressure Drop	0.04KPa@30L/min
Moisture Output	38mgH <sub>2</sub> O/L
Weight	10g
Package	200pcs/carton